




















MASTER GUIDE SCREWS ■ DRILLBITS ■ TAPS ■ INSTRUMENTS

TYPE OF SCREW	KM NO.	73-006-20	73-106-124	73-206-224	73-310-350	73-410-440	73-510-550	73-814-870	73-725-770	73-900-095	74-000-095	74-100-195
			CORTICAL				CANCELLOUS		CORTICAL	MALLEOLAR	CANCELLOUS	
diameter		1.5mm	2.0mm	2.7mm	3.5mm	3.5mm	4.0mm	4.5mm	4.5mm	6.5mm	6.5mm	6.5mm
DRILL FOR GLIDING HOLE	KM NO.	48-052	48-054	48-058	48-063	48-063	—	48-067	—	48-067	48-068	48-068
	diameter	1.5mm	2.0mm	2.7mm	3.5mm	3.5mm	NONE	4.5mm	NONE	4.5mm	*FOR SHAFT IN HARD BONE	
												
					FINE THREAD	COARSE THREAD						
DRILL FOR THREADED HOLE	KM NO.	48-051	48-052	48-054	48-057	48-054 or 48-057	48-054 or 48-057	48-061	48-061	48-061	48-061	48-061
	diameter	1.1mm	1.5mm	2.0mm	2.5mm	2.0mm or 2.5mm	2.0mm or 2.5mm	3.2mm	3.2mm	3.2mm	3.2mm	3.2mm
TAP	KM NO.	48-106	48-108	48-110	48-112	48-114	48-114	48-116	48-116	48-122	48-122	48-122
	diameter	1.5mm	2.0mm	2.7mm	3.5mm	3.5mm	3.5mm	4.5mm	4.5mm	6.5mm	6.5mm	6.5mm
					FINE THREAD	COARSE THREAD						

TYPE OF SCREW	KM NO.	73-006-20	73-106-124	73-206-224	73-310-350	73-410-440	73-510-550	73-814-870	73-725-770	73-900-095	74-000-095	74-100-195
		CORTICAL				CANCELLOUS		CORTICAL	MALLEOLAR	CANCELLOUS		
	diameter	1.5mm	2.0mm	2.7mm	3.5mm	3.5mm	4.0mm	4.5mm	4.5mm	6.5mm	6.5mm	6.5mm
COUNTERSINK	KM NO.	47-920	47-920	47-922	47-922	47-922	47-922	47-926	47-924	47-926	47-926	47-926
		diameter	1.1mm	1.1mm	2.0mm	2.0mm	2.0mm	2.0mm	4.5mm	3.2mm	4.5mm	4.5mm
DEPTH GAUGE	KM NO.	47-934	47-934	47-936	47-936	47-936	47-936	47-940	47-940	47-940	47-940	47-940
		length	24mm	24mm	50mm	50mm	50mm	50mm	100mm	100mm	100mm	100mm
QUICK COUPLING HANDLE	KM NO.	47-902	47-902	47-904	47-904	47-904	47-904	47-904	47-904	47-904	47-904	47-904
	 KM 47-904  KM 47-960			47-960	47-960	47-960	47-960	47-960	47-960	47-960	47-960	47-960
SCREWDRIVERS	KM NO.	47-950	47-950	47-953	47-953	47-953	47-953	47-954	47-954	47-954	47-954	47-954
				47-952	47-952	47-952	47-952	48-352	48-352	48-352	48-352	48-352
				48-350	48-350	48-350	48-350	48-479	48-479	48-479	48-479	48-479
				48-480	48-480	48-480	48-480					

Dimensional and material standards of instruments, bone screws and plates are according to the International Organization for Standardization (ISO).